

FALCON


ANALYTICAL SYSTEMS & TECHNOLOGY




FALCON
ANALYTICAL SYSTEMS & TECHNOLOGY



In the Lab...
In the Process



=



Data Equivalency Leads to Fluid Free Online Analyzer Calibrations

John Crandall, presenter, Joe Perron, Steve Bostic

FALCON
ANALYTICAL SYSTEMS & TECHNOLOGY



Brian Rohrback



Data Equivalency Leads to Fluid Free Analyzer Calibrations

John Crandall, presenter, Joe Perron, Steve Bostic – Falcon Analytical
Brian Rohrback - Infometrix


The age old arguments between process GC analyzer and lab GC instrument results can be mitigated by using the same kind of instrument, same method and attention to the details in both locations. Of course, the sample taken out of the process GC system is really not the same by the time it is introduced to the lab GC system. But striving to make the two units deliver “the same” result has benefits beyond settling the squabble of “who is right.”

A recent implementation in one refinery shows how to use laboratory GC as a way to calibrate the process GC using data files, not calibration fluids. This is achieved by using ASTM D7798 standard methods on both units. Expensive calibration standards and Reference Gas oil are used on the lab GC only. The resulting calibration data files from the lab GC are transferred to the process GC server. Validation occurs by use of a lab GC analyzed process sample to verify the results on the process GC. This eliminates painful operations handling cylinders of expensive calibration standards on the process GC.

Modern Processing Plant Architecture




Old School...two managers



Modern Processing Plant Architecture




New Thinking... One Manager



Lab Instruments & Process Analyzers Are Different



- Installation sites mandate design differences
 - Explosion hazard mitigation
 - Arc & spark
 - Hot surfaces
 - Instrument survival criteria in the
 - Nasty process environment vs
 - Nice laboratory setting
- All of these mandated technique differences
 - Maximum surface temperature limits (T-ratings)
 - Use of high thermal mass for temperature stability
 - Column switching vs programmed temperature operation




What if They Could Be Exactly The Same?


Data Equivalency... Regardless of Location



LAB



PROCESS





TRANSPORTABLE



Calidus 101-HT with Liquid Syringe Autosampler




- ASTM D7798 was followed
 - Split injection
 - 350°C
 - 70 nanoliter sample injected
 - Split ratio ~ 50:1
 - FID at 350°C
 - MXT-1 HT resistively heated capillary stainless steel column module
 - 320 micron ID x 0.2 micron film x 2 meter length
 - Initial temperature 40°C
 - Programmed temperature rate 2°C per second
 - Final temperature 385°C
 - Injection to injection cycle time <5 minutes




Non-injection Blank


(Autosampler - performs the injection into a waste vial and starts the PTGC cycle)



RTMIX (Autosampler - injects the retention time standard mix, C₅ to C₄₄)



RGO (Autosampler - injects the Reference Gas Oil sample)



RGO SimDist Result using the Autosampler




All boiling points are within the tolerance of ASTM D-7798.

RGO (deg. F)				
<u>% OFF</u>	<u>Accepted</u>	<u>Measured</u>	<u>Difference</u>	<u>D-7798 Tolerance</u>
IBP	239	236.4	-2.6	12.6
10	349	346.1	-2.9	8
20	435	434.3	-0.7	9
30	499	498.3	-0.7	8.6
40	552	552.2	0.2	7.7
50	594	593.8	-0.2	7.7
60	629	628.7	-0.3	7.7
70	669	669.3	0.3	7.7
80	712	713.6	1.6	7.7
90	764	767.1	3.1	7.7
FBP	887	897.9	10.9	21.2

Samples Run with Autosampler - LBD, MBD, HBD, Kerosene



The runs demonstrate the similarity and, of course, the differences in the four samples.




FALCON
ANALYTICAL SYSTEMS & TECHNOLOGY



Process Use




- Injections Using a Valco™ Liquid Injection Sample Valve
- Boiling Point Analysis using Simulated Distillation



Fast GC with Rotary Liquid Injection Valve




- ASTM D7798 was followed
 - Sample inlet via heated Valco rotary valve
 - Maintained at 225°C
 - Syringe loaded with capillary outlet to maintain pressure in the loop
 - 3 second delay was employed for pressure balance reasons (now proven unnecessary)
 - Split injection
 - 350°C
 - 60 nanoliter sample injected
 - Split ratio ~ 50:1
 - FID at 350°C
 - MXT-1 HT resistively heated capillary stainless steel column module
 - 320 micron ID x 0.2 micron film x 2 meter length
 - Initial temperature 40°C
 - Programmed temperature rate 2°C per second
 - Final temperature 385°C
 - Injection to injection cycle time <5 minutes




Non-injection Blank


(Valve - actuator is turned off and the PTGC cycle is started)



RTMIX (Valve - cleaned and loaded using a syringe with enough backpressure to keep components from flashing too soon)



RGO (Valve- injects the Reference Gas Oil sample)




RGO SimDist Results (Valve)



All boiling points are within the tolerance of ASTM D-7798.

RGO (deg. F)				
<u>% OFF</u>	<u>Accepted Value</u>	<u>Measured Value</u>	<u>Difference</u>	<u>D-7798 Tolerance</u>
IBP	239	237.2	-1.8	12.6
10	349	347.1	-1.9	8
20	435	437.1	2.1	9
30	499	496.9	-2.1	8.6
40	552	556.5	4.5	7.7
50	594	596.6	2.6	7.7
60	629	629.9	0.9	7.7
70	669	671.0	2.0	7.7
80	712	714.5	2.5	7.7
90	764	767.0	3.0	7.7
FBP	887	889.0	2.0	21.2

Samples Run with Valve - LBD, MBD, HBD, Kerosene





Lab and Process Comparison




Autosampler Inlet



=

Rotary Valve Inlet



RGO (Autosampler Minus Valve Results)



RGO Autosampler		RGO Valve		
<u>% OFF</u>	<u>BP(F)</u>	<u>% OFF</u>	<u>BP(F)</u>	<u>Difference</u>
IBP	236.4	IBP	237.2	0.8
10	346.1	10	347.1	1.0
20	434.3	20	437.1	2.8
30	498.3	30	496.9	-1.4
40	552.2	40	556.5	4.3
50	593.8	50	596.6	2.8
60	628.7	60	629.9	1.2
70	669.3	70	671	1.7
80	713.6	80	714.5	0.9
90	767.1	90	767	-0.1
FBP	897.9	FBP	889	-8.9

RGO Autosampler Minus Valve (°F)



% OFF	LBD		MBD		HBD		Kero	
	Allowable		Allowable		Allowable		Allowable	
	AVG	Delta	AVG	Delta	AVG	Delta	AVG	Delta
IBP	-0.4	11.8	0.6	11.8	-2.0	11.8	-2.3	11.8
5	0.2	7.2	0.2	7.2	-1.4	7.2	-0.2	7.2
10	-1.5	7.7	0.0	7.7	0.8	7.7	0.8	7.7
20	-1.2	8.5	-0.4	8.5	1.9	8.5	1.9	8.5
30	3.9	7.8	4.8	7.8	5.4	7.8	2.2	7.8
40	7.3	7.7	2.6	7.7	4.1	7.7	2.1	7.7
50	4.9	7.7	-0.5	7.7	-3.2	7.7	2.8	7.7
60	-3.0	7.7	-1.7	7.7	-2.2	7.7	2.6	7.7
70	-2.3	7.7	-3.4	7.7	-0.9	7.7	4.4	7.7
80	-1.3	7.7	-3.7	7.7	-2.5	7.7	6.9	7.7
90	-1.9	7.7	-2.9	7.7	-3.7	7.7	9.1	7.7
95	-3.9	9.0	-3.1	9.0	-2.5	9.0	7.1	9.0
FBP	-13.8	21.2	-12.5	21.2	-13.1	21.2	-16.5	21.2

Note:

8 replicates of each sample type were run with each injection method. These results were averaged to obtain the data in this table. Allowable delta is determined according to ASTM D2887 reproducibility rules

4 Samples, 8 Replicate Statistics (°F) (Autosampler and Valve)



% OFF	Autosampler								Valve							
	LBD		MBD		HBD		Kero		LBD		MBD		HBD		Kero	
	AVG	STDEV	AVG	STDEV	AVG	STDEV	AVG	STDEV	AVG	STDEV	AVG	STDEV	AVG	STDEV	AVG	STDEV
IBP	228.0	0.5	225.7	0.2	238.0	0.2	224.8	0.3	228.3	2.4	225.1	3.5	239.9	0.3	227.1	0.5
5	329.1	0.1	329.5	0.1	339.9	0.3	292.5	0.1	329.0	0.2	329.3	0.7	341.4	0.4	292.7	0.3
10	363.3	0.2	370.1	0.1	384.3	0.2	319.4	0.1	364.8	0.6	370.1	1.0	383.6	0.1	318.6	0.2
15	391.1	0.2	404.0	0.1	415.6	0.4	333.0	0.1	391.9	0.3	404.9	1.4	416.9	0.5	331.8	0.3
20	416.3	0.2	432.5	0.1	442.4	0.4	345.5	0.0	417.5	0.2	432.9	1.5	440.6	0.3	343.6	0.2
25	435.7	0.2	456.6	0.1	462.8	0.5	353.6	0.1	435.9	0.2	453.2	0.7	461.0	0.6	351.8	0.4
30	455.4	0.2	480.2	0.1	485.0	0.6	362.8	0.1	451.4	0.1	475.5	0.8	479.7	0.3	360.6	0.3
35	473.5	0.2	495.1	0.1	503.3	0.7	373.3	0.1	469.7	0.1	490.5	1.5	498.6	0.5	370.9	0.3
40	489.1	0.2	514.2	0.1	520.5	0.7	384.4	0.1	481.9	0.1	511.6	1.1	516.4	0.5	382.2	0.5
45	505.4	0.2	529.8	0.1	538.1	0.8	389.3	0.1	500.4	0.3	529.6	1.0	538.8	0.7	386.4	0.4
50	520.7	0.2	548.2	0.0	555.7	1.0	398.8	0.1	515.8	0.3	548.6	0.5	559.0	0.6	396.0	0.5
55	536.7	0.2	565.0	0.1	575.0	1.0	409.2	0.1	536.6	0.4	567.0	0.9	575.1	0.4	406.5	0.5
60	551.9	0.2	579.7	0.1	589.7	1.1	419.0	0.1	554.9	0.7	581.4	0.9	591.8	0.6	416.4	0.6
65	571.1	0.3	598.9	0.1	605.7	1.3	425.7	0.1	573.2	0.4	600.7	0.3	608.6	0.7	422.2	0.5
70	585.4	0.3	616.0	0.2	626.3	1.3	436.2	0.1	587.7	0.6	619.4	0.9	627.2	0.5	431.8	0.7
75	603.0	0.2	636.1	0.2	647.2	1.4	447.6	0.1	604.3	0.7	639.8	0.9	649.4	0.5	442.1	0.6
80	624.4	0.3	658.0	0.2	669.2	1.5	458.4	0.1	625.8	0.5	661.7	1.0	671.7	0.6	451.5	0.4
85	647.5	0.3	682.7	0.2	694.7	1.6	471.2	0.1	649.6	0.7	686.2	1.2	696.1	0.7	463.9	1.3
90	673.9	0.2	712.3	0.2	726.4	1.7	488.9	0.1	675.8	1.7	715.2	1.0	730.1	1.2	479.8	0.7
95	711.0	0.3	751.9	0.3	773.3	1.9	511.3	0.1	714.9	1.9	755.0	1.9	775.9	2.0	504.2	3.2
FBP	793.1	1.3	828.6	1.1	856.2	2.3	576.7	0.2	806.9	8.9	841.1	12.7	869.3	12.3	593.2	26.8
Avg. St. Dev	0.3		0.2		1.0		0.1		1.0		1.7		1.2		1.9	

4 Sample, All 16 Replicates Combined Statistics (°F)

(Autosampler and Valve)



Combined									
	LBD		MBD		HBD		Kero		
<u>% OFF</u>	<u>AVG</u>	<u>STDEV</u>	<u>AVG</u>	<u>STDEV</u>	<u>AVG</u>	<u>STDEV</u>	<u>AVG</u>	<u>STDEV</u>	<u>STDEV</u>
IBP	228.1	1.7	225.4	2.4	238.9	1.0	225.9	1.3	
5	329.0	0.2	329.4	0.5	340.7	0.8	292.6	0.3	
10	364.1	0.9	370.1	0.7	384.0	0.4	319.0	0.4	
15	391.5	0.5	404.4	1.1	416.2	0.8	332.4	0.6	
20	416.9	0.7	432.7	1.1	441.5	1.0	344.6	1.0	
25	435.8	0.2	454.9	1.8	461.9	1.1	352.7	1.0	
30	453.4	2.0	477.8	2.5	482.4	2.8	361.7	1.2	
35	471.6	1.9	492.8	2.6	500.9	2.5	372.1	1.3	
40	485.5	3.8	512.9	1.5	518.4	2.2	383.3	1.2	
45	502.9	2.6	529.7	0.7	538.4	0.8	387.8	1.5	
50	518.3	2.6	548.4	0.4	557.3	1.9	397.4	1.5	
55	536.7	0.3	566.0	1.2	575.1	0.7	407.8	1.5	
60	553.4	1.6	580.5	1.1	590.8	1.4	417.7	1.4	
65	572.1	1.1	599.8	1.0	607.1	1.8	424.0	1.9	
70	586.6	1.3	617.7	1.8	626.8	1.1	434.0	2.3	
75	603.6	0.8	638.0	2.0	648.3	1.6	444.9	2.9	
80	625.1	0.8	659.8	2.1	670.4	1.7	454.9	3.6	
85	648.6	1.2	684.4	2.0	695.4	1.4	467.5	3.9	
90	674.9	1.5	713.7	1.6	728.3	2.4	484.3	4.7	
95	712.9	2.4	753.5	2.1	774.6	2.3	507.8	4.3	
FBP	800.0	9.4	834.8	10.8	862.7	10.9	585.0	20.2	
Avg. Std. Dev.	1.8		1.9		1.9		2.7		

FALCON
ANALYTICAL SYSTEMS & TECHNOLOGY




LineUp™

Retention Time Alignment

and

Pirouette®

Chemometrics



Autosampler and Valve Comparison

LineUp and Pirouette Were Used for the Correction of Retention Time Differences and Area Normalization







FALCON
ANALYTICAL SYSTEMS & TECHNOLOGY



**Data
Equivalency**

All 16 Autosampler and Valve Results Overlaid after Alignment and Normalization




PCA Scores Plot – 16 Runs Per Sample Type

(8 Autosampler and 8 Valve, all results within ASTM D-7798 tolerances)



Y-axis amplified x50







- Autosampler
- Sample Valve

SimDist D7798: lab versus on-line



Laboratory




- 1. Run a Blank**
to insure a clean system
- 2. Run the n-Paraffin Standard**
relates retention time to temperature
- 3. Run the Reference**
to check results with a known
- 4. Run Samples**



SimDist D7798: lab versus on-line



On-Line




Pass



1. Run a Blank
to insure a clean system
2. Use the **same n-Paraffin Standard**
from the run in the laboratory
3. Use the **RGO or Process Sample**
as an alignment target
4. Run Samples

- The exact same GC model and methods can be used in the lab and the process
- What does this mean?
 - Lab GC data can completely ***eliminate*** the need for time consuming, complex work and ***costly calibration standards*** used on process GC's
 - The normal calibration process is: Blank → RTMIX → RGO
 - With chromatographic alignment, the process GC can be calibrated with ONLY ONE BLANK RUN using lab GC data files




Data Equivalency...

Regardless of Location the Result Are the Same



LAB



PROCESS



TRANSPORTABLE



Questions?



*Thank you for your
interest & attention!*

